
MultiMotion

Dynamic Corrective Joint

Instructions for use

Basko
Healthcare

Multi**Motion**

NL	- Technische gebruiksaanwijzing	5
DE	- Technische Gebrauchsanleitung.....	19
ENG	- Technical manual.....	33
FR	- Manuel technique.....	47
ES	- Instrucciones técnicas.....	61
IT	- Istruzioni para el usuario	75
SI	- Tehnična navodila.....	89
CZ	- Technické pokyny.....	103
SE	- Tekniska anvisningar.....	117
DK	- Tekniske instruktioner	131
NO	- Tekniske instruksjoner.....	145
FI	- Tekniset ohjeet.....	159

Multi**Motion**

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Multi**Motion**

Please read these user instructions carefully before starting to use the product. Pay particular attention that you follow the safety instructions.

Usage Guidelines

- The product must be fitted by or under the supervision of a competent professional.
- Special precautions must always be taken for patients with a visual impairment, cognitive limitations and/or reduced sensibility in the upper and/or lower extremity.
- Changes, modifications and adjustments to the product that are not described in these instructions, are done under the responsibility of the person carrying them out.

Intended use of the product

Orthotics with a MultiMotion dynamic correction system joint - hereinafter called the 'system joint' - are used for the treatment of dynamic corrective limb contractures in accordance with the Low-Load Prolonged Stretch (LLPS) method. This involves applying a low force to a limb over a longer period of time in order to correct or prevent the contracture. The system joint is to be used exclusively in orthotics for the upper or lower extremities and should not be weight bearing.

Indications / Contraindications

The MultiMotion system joints can be used with children and adults for the treatment of dynamically correctable contractures and/or for their prevention.

The origin of contracture is a multitude of indications that can be given by a doctor, among others:

Neurological indications:

- Apoplexy (CVA)
- Cerebral paresis
- Dystrophy
- Multiple sclerosis
- Paresis
- Myelomeningocele
- Skull and brain trauma

Orthopaedic indications:

- Amputations
- Band rupture
- Fracture
- Post-op after placing endoprostheses, etc.

Contraindications

Structural contractures:

- Abnormal arthrokinematics
- Ankylosis
- Arthrodesis
- Fibrosis
- Ossification

General safety instructions

Risk of injury as a result of incorrect use; the professional should verify correct functioning of the system joints prior to delivery of the orthotics to the user. The patient and his/her carer(s) should also be instructed how to put on the orthotics and have the joint operation explained to them. It should also be made clear to them that the spring tension setting should never be changed by anybody else than the professional or therapist, unless the treating physician decides otherwise. The joints should not come into contact with water.

Please note: the operating principle is not the higher the force the faster the R.O.M. gain. On the contrary, excessive force can result in retraction from muscle tissue or micro lesions in tissue.

Comment

Specialists applying and using the MultiMotion system joint are expected to use and apply their professional judgement when selecting, fitting, adjusting and adequately educating the patient or carer, in order to minimise the potential risks to each individual patient.

These risks may include the mentioned contraindications, but may also be associated with the unique patient characteristics or circumstances of the patient's carer. In addition, there are possible risks of creating micro lesions, retractions in the muscle belly, inducing spasms, pressure spots or other.

Materials used

The MultiMotion system joints incorporate components made of various types of materials, including:

- Steel/stainless steel
- Plastic

Maintenance recommendation

For optimum operation, keep the tension spring free from dirt and other external influences that could affect its operation. We recommend checking the joints every 6 months for functionality and possible wear and tear.

Cleaning instructions




The joint parts of the system joint can be cleaned but are not resistant to aggressive agents such as acetone or similar. Try to keep the joint free from dust and dirt as much as possible.

Package contents

- MultiMotion system joint
- Allen key
- Torx Key
- MultiMotion technical instructions
- MultiMotion patient documentation form

The system joint is available in several variants. Check the package label to make sure the shown article number corresponds with the contents of the package.

Components



1. Bar adapter with screws
2. Point of deflection
3. Central housing
4. Central housing screw
5. Force output indicator


6. Force output indicator disc
7. Spring tension adjustment screw
8. Lock / unlock lever
9. Flexion / extension stop screw

Dimensions MultiMotion Medium



MODULAR (#MM04.2M-12)


MODULAR T-BAR (#MM04.2MT-12)

INTEGRATED (#MM04.2B-12)

INTEGRATED T-BAR (#MM04.2BT-12)

Dimensions MultiMotion Large

MODULAR (#MM05.2M-16)

MODULAR T-BAR (#MM05.2MT-16)


INTEGRATED (#MM05.2B-16)

INTEGRATED T-BAR (#MM05.2BT-16)

Product & parts MultiMotion Medium


PRODUCT

NO.	ART. NO.	DESCRIPTION	NUMBER OF
A1	MM04.2M-12	MULTIMOTION MODULAR MEDIUM	UNIT
A2	MM04.2B-12	MULTIMOTION INTEGRATED MEDIUM	UNIT
A3	MM04.2MT-12	MULTIMOTION MODULAR T-BAR MEDIUM	UNIT
A4	MM04.2BT-12	MULTIMOTION INTEGRATED T-BAR MEDIUM	UNIT
A5	MM00-01	ALLEN KEY - 2.5 MM	UNIT
A6	MM00-01T	TORX KEY - T8	UNIT

PARTS (SEE COVER)

NO.	ART. NO.	DESCRIPTION	NUMBER OF
B	MM04.2-10	CENTRAL HOUSING SCREW	UNIT
C	MM04.2-20	HOUSING	UNIT
D	MM04.2-30	WASHER	UNIT
E	MM04.2-40	FORCE OUTPUT INDICATOR DISC	UNIT
F	MM04.2-50	SPRING CASSETTE	UNIT
G	MM04.2-60	SLIDE RING	UNIT
H	MM04.2-70	FLEXION / EXTENSION STOPS	SET (2X)
I	MM00-02	SCREWS (PH-2) FOR SYSTEM BARS	SET (4X)

Product & parts MultiMotion Large


PRODUCT

NO.	ART. NO.	DESCRIPTION	NUMBER OF
A1	MM05.2M-16	MULTIMOTION MODULAR LARGE	UNIT
A2	MM05.2B-16	MULTIMOTION INTEGRATED LARGE	UNIT
A3	MM05.2MT-16	MULTIMOTION MODULAR T-BAR LARGE	UNIT
A4	MM05.2BT-16	MULTIMOTION INTEGRATED T-BAR LARGE	UNIT
A5	MM00-01	ALLEN KEY - 2.5 MM	UNIT
A6	MM00-01T	TORX KEY - T8	UNIT

PARTS (SEE COVER)

NO.	ART. NO.	DESCRIPTION	NUMBER OF
B	MM05.2-10	CENTRAL HOUSING SCREW	UNIT
C	MM05.2-20	HOUSING	UNIT
D	MM05.2-30	WASHER	UNIT
E	MM05.2-40	FORCE OUTPUT INDICATOR DISC	UNIT
F	MM05.2-50	SPRING CASSETTE	UNIT
G	MM05.2-60	SLIDE RING	UNIT
H	MM05.2-70	FLEXION / EXTENSION STOPS	SET (2X)
I	MM00-02	SCREWS (PH-2) FOR SYSTEM BARS	SET (4X)

Overview system bars & accessories

NO.	ART. NO.	DESCRIPTION	NUMBER OF
I	BAR-S-12	SYSTEM BAR, STRAIGHT 200 X 12 X 3.2 MM	2 PIECES
	BAR-S-16	SYSTEM BAR, STRAIGHT 200 X 16 X 4 MM	2 PIECES
II	BAR-O-12	SYSTEM BAR, CURVED 350 X 12 X 3.2 MM	2 PIECES
	BAR-O-16	SYSTEM BAR, CURVED 450 X 16 X 4 MM	2 PIECES
III	BAR-T-12	SYSTEM BAR, 150 X 12 X 3.2 MM	2 PIECES
	BAR-T-16	SYSTEM BAR, 150 X 16 X 3.2 MM	2 PIECES
IV	AK-100	BASIC ALIGNMENT T-BAR	UNIT
	AK-102	MULTIMOTION ALIGNMENT UNIT	SET

Joint selection

The MultiMotion system joint is available in two sizes, 'Medium' and 'Large' and are both suitable for children, adolescents and adults. We generally recommend the medium joint for elbow and wrist joints and the large joint for knee and ankle joints, but this depends on the morphology and condition of the patient.

MultiMotion is universally applicable in the treatment of flexion or extension contractures. The removable spring cassette allows the gripping point and spring force to be adjusted in the range of movement and read off in any angular position.

Both system joints have a total range of movement of -90° to 90° with an adjustable flexion or extension stop. The spring tension is infinitely adjustable with the supplied Allen key and the joint itself can be locked and unlocked in almost any position.

Please note: ensure that the system joint is positioned correctly, i.e. the MultiMotion logo must be legible on the housing (*not upside down*).

The joints can be used at both the left and right side and can be set with force directed in flexion or extension. This direction of force depends on the patient's body side and the specific joint on which the joint is applied. The spring cassette is painted red on one side and yellow on the other. The chart below shows the force direction per joint, body side and spring cassette colour. As standard, the joint is supplied with the red side of the spring cassette facing up.

Adjusting flexion/extension

BODY JOINT	BODY SIDE	FORCE IN FLEXION (PLANTAR FLEXION)		FORCE IN EXTENSION (DORSIFLEXION)	
		YELLOW	RED	YELLOW	RED
WRIST	LEFT	X			X
	RIGHT		X	X	
ELBOW	LEFT		X	X	
	RIGHT	X			X
KNEE	LEFT	X			X
	RIGHT		X	X	
ANKLE	LEFT	X			X
	RIGHT		X	X	

Changing the direction of force (see cover)

The description below assumes that the system joint has been delivered factory-delivered and therefore the spring cassette is mounted with the red side facing up.

Step 1

Loosen the central screw on the housing using the torx wrench provided. Now remove the screw (*Fig. 1*).

Step 2

Next, lift up the metal housing over the central shaft. Now a white plastic washer and the force indicator disc with the red printing are visible (*Fig. 2*). Lift both up from the spring cassette and put them safely aside.

Step 3

Now lift the spring cassette up over the central axis from the joint and turn it over completely so the yellow cassette side becomes visible (*Fig. 3*).

Step 4

Then place the cassette back over the central axis but with the **yellow** side up. Essential here is the 'smiley' (☺) marking on the cassette is positioned at '6 o'clock' (*Fig. 4*).

Step 5

Next, place the force indicator disc on the spring cassette - the yellow print facing up (!) - followed by the white plastic washer. Make sure the digit 0 on the force indicator disc is aligned with Force output indicator marking (⑥) on the joint (*Fig. 5*)!

Step 6

Next, reposition the housing and secure it with the center screw (*Fig. 6*).

Mounting the joint on system bars


To mount the bars on the modular system joint, place the bars - with the countersinks facing down - over the bushings of the joint. Fasten the bars with the supplied screws and a PH-2 screwdriver and secure them with Loctite 243.

Aligning the joint

In order for the joints to function correctly, it is important that the system joint and the carriage are aligned correctly one to another, and are mounted on an orthotics construction with sufficient torsion rigidity. During orthotic production, use the Basic Alignment T-Bar (AK-100) + MultiMotion alignment unit (AK-102).

Bending the joint

Proximal and distal to the joint are small notches (②). The joint bars must be bent horizontally between these bending points. Place the joint horizontally with mounted bar(s) between the bench screw. Do not use bending irons but push/pull the joint to the desired angle so that the joint bends on the line of the notch.

Calibrate joint

This should be done on the plaster positive. To do this, first remove the spring cassette (see steps 1 to 3 'Changing the direction of force') and the flexion / extension stops at the rear. Then position the orthotics on the plaster positive and push the locking / unlocking lever to the horizontal position to lock the angle adjustment of the orthotics. Then remove the orthotics from the cast and fit the flexion / extension stops, adjusted to the set angle.

Please note: Check that the visible color on the spring cassette is in accordance with the intended flexion/extension force direction. Place the force output indicator disc with the same color as the spring cassette this, but such that the digit 0 is aligned with the proximal force output indicator marking (6) on the brake housing of the system joint. After this, insert the washer and housing and secure it with the central housing screw, then unlock the hinge.

Locking the joint

In order to more easily position and remove the orthotics, the joint can be locked at virtually each desired angle. The release and locking mechanism is located to the side of the lower joint half. In order to lock the joint, the catch should be pushed forward/up.

Setting flexion and extension stops

To ensure correct operation, it is mandatory to set the flexion and extension stops. Failure to do so may result in excessive strain on the joint. There are two screws at the back of the functional MultiMotion joint to set the range of motion. The standard stop setting in both functional joints is 0° flexion and 90° extension. With the 'Medium' model, the stops can be moved in 15° increments and with the 'Large' model in 12° increments. Use the supplied roundhead screwdriver to loosen the flexion and extension screws, move them as needed and then retighten them.

CAUTION

BEWARE OF RISKS OF INJURY BEFORE UNLOCKING THE JOINT. THE SET SPRING TENSION MUST BE TAKEN INTO ACCOUNT. IF THIS IS HIGH, IT IS ADVISABLE TO HAVE A CARER ASSIST WITH UNLOCKING IN ORDER TO GRADUALLY GUIDE THE ORTHOTICS TOWARDS THE SET FLEXION OR EXTENSION FORCE.

Setting the spring tension

On both sides of the system joint, there is an Allen wrench with which the spring tension can be adjusted.

Please note: there should be minimum tension (0) during the fitting phase. It should only be adjusted during rehabilitation, taking into account that the patient may have to get used to wearing the orthotics first. The duration of habituation is determined by the sensibility and is therefore patient-dependent and can vary from hours to days and is determined by the treating specialist/therapist.

Patient documentation form

A patient documentation form is included in the package on which the progression of the active and passive range of motion, spring tension and wearing time can be recorded. Be sure to give this form to the patient/carer or therapist concerned.

Liability

Basko Healthcare as the manufacturer is solely liable if the product is used in accordance with the prescribed conditions and the purpose for which it is intended. Basko Healthcare recommends handling the product in accordance with applicable rules and performing maintenance as described in these user instructions. Non-compliance with the above may cause the MultiMotion joints to not function or function inadequately and may impact the guarantee. Changes, adjustments and modifications to the product not described in these instructions, are the responsibility and liability of the person who performs them.

CE conformity

These products meet the requirements in Directive (EU) 2017/745 for medical devices. These products have been included in Class I based on classification criteria for medical devices as per annex VIII of the Directive. Basko Healthcare has therefore provided the EU Declaration of Conformity under its own sole responsibility in accordance with Annex XI of the Directive.

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