BASKO HEALTHCARE
 Pieter Lief tinckweg 16
 NL-1505 HX, Zaandam
 

Basko Health care

Pieter Lief tinckweg 16 | NL-1505 HX Zaandam
 Tel: +31 (0) 75 - 613 15 13 | Fax: +31 (0) 75 - 612 63 73
 E-mail: sales@basko.com | Internet: www.basko.com

basko.com

MM-9950 / Stand: 04/2022

MultiMotion

Dynamic Corrective Joint

Instructions for use

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MultiMotion

Please read these user instructions carefully before starting to use the product. Pay particular attention that you follow the safety instructions.

Usage Guidelines

- The product must be fitted by or under the supervision of a competent professional.
- Special precautions must always be taken for patients with a visual impairment, cognitive limitations and/or reduced sensibility in the upper and/or lower extremity.
- Changes, modifications and adjustments to the product that are not described in these instructions, are done under the responsibility of the person carrying them out.

Intended use of the product

Orthotics with a MultiMotion dynamic correction system joint - hereinafter called the 'system joint' - are used for the treatment of dynamic corrective limb contractures in accordance with the Low-Load Prolonged Stretch (LLPS) method. This involves applying a low force to a limb over a longer period of time in order to correct or prevent the contracture. The system joint is to be used exclusively in orthotics for the upper or lower extremities and should not be weight bearing.

Indications / Contraindications

The MultiMotion system joints can be used with children and adults for the treatment of dynamically correctable contractures and/or for their prevention.

The origin of contracture is a multitude of indications that can be given by a doctor, among others:

Neurological indications:

- Apoplexy (CVA)
- Cerebral paresis
- Dystrophy
- Multiple sclerosis
- Paresis
- Myelomeningocele
- Skull and brain trauma

Orthopaedic indications:

- Amputations
- Band rupture
- Fracture
- Post-op after placing endoprostheses, etc.

Contraindications

Structural contractures:

- Abnormal arthrokinematics
- Ankylosis
- Arthrodesis
- Fibrosis
- Ossification

General safety instructions

Risk of injury as a result of incorrect use; the professional should verify correct functioning of the system joints prior to delivery of the orthotics to the user. The patient and his/her carer(s) should also be instructed how to put on the orthotics and have the joint operation explained to them. It should also be made clear to them that the spring tension setting should never be changed by anybody else than the professional or therapist, unless the treating physician decides otherwise. The joints should not come into contact with water. Furthermore, it is important that the set force is set as low as possible without affecting the corrective force.

Please note: the operating principle is not the higher the force the faster the R.O.M. gain. On the contrary, excessive force can result in retraction from muscle tissue or micro lesions in tissue.

Comment

Specialists applying and using the MultiMotion system joint are expected to use and apply their professional judgement when selecting, fitting, adjusting and adequately educating the patient or carer, in order to minimise the potential risks to each individual patient.

These risks may include the mentioned contraindications, but may also be associated with the unique patient characteristics or circumstances of the patient's carer. In addition, there are possible risks of creating micro lesions, retractions in the muscle belly, inducing spasms, pressure spots or other.

Materials used

The MultiMotion system joints incorporate components made of various types of materials, including:

- Steel/stainless steel
- Plastic

Maintenance recommendation

The tension spring in the system joint is already greased with PTFE spray during assembly. For optimum operation, keep the tension spring free from dirt and other external influences that could affect its operation. We recommend checking the joints every 6 months for functionality and possible wear and tear.

Cleaning instructions




The joint parts of the system joint can be cleaned but are not resistant to aggressive agents such as acetone or similar. Try to keep the joint free from dust and dirt as much as possible.

Package contents

- MultiMotion system joint
- Allen key
- MultiMotion technical instructions
- MultiMotion patient documentation form

The system joint is available in several variants. Check the package label to make sure the shown article number corresponds with the contents of the package.

Components

- | | |
|----------------------------|------------------------------------|
| 1. Bar adapter with screws | 6. Force output indicator disc |
| 2. Point of deflection | 7. Spring tension adjustment screw |
| 3. Central housing | 8. Lock / unlock lever |
| 4. Central housing screw | 9. Flexion / extension stop screw |
| 5. Force output indicator | 10. Flexion / extension settings |

Dimensions MultiMotion Medium

MODULAR (#MM04M-12)


MODULAR T-BAR (#MM04MT-12)

INTEGRATED (#MM04B-12)

INTEGRATED T-BAR (#MM04BT-12)


Dimensions MultiMotion Large

MODULAR (#MM05M-16)

MODULAR T-BAR (#MM05MT-16)


INTEGRATED (#MM05B-16)

INTEGRATED T-BAR (#MM05BT-16)

Product & parts MultiMotion Medium


NO.	ART. NO.	DESCRIPTION	NUMBER OF
A1	MM04M-12	MULTIMOTION MODULAR MEDIUM	UNIT
A2	MM04B-12	MULTIMOTION INTEGRATED MEDIUM	UNIT
A3	MM04MT-12	MULTIMOTION MODULAR T-BAR MEDIUM	UNIT
A4	MM04BT-12	MULTIMOTION INTEGRATED T-BAR MEDIUM	UNIT
A5	MM00-01	SOCKET SPANNER - 2.5 MM	UNIT

SEE COVER

B	MM04-10	CENTRAL HOUSING SCREW	UNIT
C	MM04-20	HOUSING	UNIT
D	MM04-30	WASHER	UNIT
E	MM04-40	FORCE OUTPUT INDICATOR DISC	UNIT
F1	MM04-50	COVER TENSION SPRING-LOADED CASSETTE	UNIT
F2		TENSION SPRING-LOADED CASSETTE 4.1 NM	
G	MM04-60	SLIDE RING	UNIT
H	MM04-70	FLEXION / EXTENSION STOPS	SET (2X)
I	MM00-02	SCREWS (PH-2) FOR SYSTEM BARS	SET (4X)

Product & parts MultiMotion Large


NO.	ART. NO.	DESCRIPTION	NUMBER OF
A1	MM05M-16	MULTIMOTION MODULAR LARGE	UNIT
A2	MM05B-16	MULTIMOTION INTEGRATED LARGE	UNIT
A3	MM05MT-16	MULTIMOTION MODULAR T-BAR LARGE	UNIT
A4	MM05BT-16	MULTIMOTION INTEGRATED T-BAR LARGE	UNIT
A5	MM00-01	SOCKET SPANNER - 2.5 MM	UNIT

SEE COVER

B	MM05-10	CENTRAL HOUSING SCREW	UNIT
C	MM05-20	HOUSING	UNIT
D	MM05-30	WASHER	UNIT
E	MM05-40	FORCE OUTPUT INDICATOR DISC	UNIT
F1	MM05-50	COVER TENSION SPRING-LOADED CASSETTE	UNIT
F2		TENSION SPRING-LOADED CASSETTE 8.4 NM	
G	MM05-60	SLIDE RING	UNIT
H	MM05-70	FLEXION / EXTENSION STOPS	SET (2X)
I	MM00-02	SCREWS (PH-2) FOR SYSTEM BARS	SET (4X)

Overview system bars & accessories

NO.	ART. NO.	DESCRIPTION	NUMBER OF
I	BAR-S-12	SYSTEM BAR, STRAIGHT 200 X 12 X 3.2 MM	2 PIECES
	BAR-S-16	SYSTEM BAR, STRAIGHT 200 X 16 X 4 MM	2 PIECES
II	BAR-O-12	SYSTEM BAR, CURVED 350 X 12 X 3.2 MM	2 PIECES
	BAR-O-16	SYSTEM BAR, CURVED 450 X 16 X 4 MM	2 PIECES
III	BAR-T-12	SYSTEM BAR, 150 X 12 X 3.2 MM	2 PIECES
	BAR-T-16	SYSTEM BAR, 150 X 16 X 3.2 MM	2 PIECES
IV	AK-100	BASIC ALIGNMENT T-BAR	UNIT
	AK-102	MULTIMOTION ALIGNMENT UNIT	SET

Joint selection

The MultiMotion system joint is available in two sizes, 'Medium' and 'Large' and are both suitable for children, adolescents and adults. We generally recommend the medium joint for elbow and wrist joints and the large joint for knee and ankle joints, but this depends on the morphology and condition of the patient.

MultiMotion is universally applicable in the treatment of flexion or extension contractures. The removable spring cassette allows the gripping point and spring force to be adjusted in the range of movement and read off in any angular position.

Both system joints have a total range of movement of -95° to 95° with an adjustable flexion or extension stop. The spring tension is infinitely adjustable with the supplied Allen key and the joint can be locked and unlocked in almost any position.

Please note: ensure that the system joint is positioned correctly, i.e. the MultiMotion logo must be legible on the housing (*not upside down*).

The joints can be used at both the left and right side and can be set with force in flexion or extension direction. This direction of force depends on the side of the body and the patient's joint for which the system joint is used. The tension spring is painted red on one side and yellow on the other. The table below shows the force direction per joint, body side and spring colour. As standard, the joint is delivered with the red side of the tension spring facing upwards in the cassette.

Adjusting flexion/extension

BODY JOINT	BODY SIDE	FORCE IN FLEXION (PLANTAR FLEXION)		FORCE IN EXTENSION (DORSIFLEXION)	
		YELLOW	RED	YELLOW	RED
WRIST	LEFT	X			X
	RIGHT		X	X	
ELBOW	LEFT		X	X	
	RIGHT	X			X
KNEE	LEFT	X			X
	RIGHT		X	X	
ANKLE	LEFT	X			X
	RIGHT		X	X	

Changing the direction of force (See right cover)

The following description assumes that the system joint is supplied new and therefore the tension spring is mounted with the red side facing upwards.

Step 1

Loosen the central screw on the housing using the roundhead screwdriver. Now remove the screw (Fig. 1).

Step 2

Lift up the housing over the central axis and lay it aside. Now a white washer and the force indicator disc with the red printing are visible (Fig. 2). Lift both also from the spring cassette and put them safely aside. The spring-loaded cassette including the spring can now be moved up from the axis over the central joint axis.

Step 3

Use a flat screwdriver to remove the cover plate of the spring cassette and put it aside (Fig. 3). The red side of the tension spring is now visible (Fig. 4).

Step 4

Lift the cassette up from the joint and turn it over completely so that the yellow spring side is visible. Then place the cassette back over the central axis on the joint (Fig. 5). **Please note:** Make sure that the slide ring under the spring cassette is correctly positioned before replacing it.

Step 5

Replace the cover plate on the spring cassette (Fig. 6) and place the force indicator disc, yellow print side up, on the cover plate. (Fig. 7). Align the 0 on the indicator disc with the force output indicator mark (6) on the brake housing.

Step 6

Replace the washer and the cap correctly and screw them back on (Fig. 8).

Mounting the joint on system bars


To mount the bars on the modular system joint, place the bars - with the countersinks facing down - over the bushings of the joint. Fasten the bars with the supplied screws and a PH-2 screwdriver and secure them with Loctite 243.

Aligning the joint

In order for the joints to function correctly, it is important that the system joint and the carriage are aligned correctly one to another, and are mounted on an orthotics construction with sufficient torsion rigidity. During orthotic production, use the Basic Alignment T-Bar (AK-100) + MultiMotion alignment unit (AK-102).

Bending the joint

Proximal and distal to the joint are small notches (②). The joint bars must be bent horizontally between these bending points. Place the joint horizontally with mounted bar(s) between the bench screw. **Do not** use bending irons but push/pull the joint to the desired angle so that the joint bends on the line of the notch.

Calibrate joint

This should be done on the plaster positive. To do this, first remove the spring cassette (see steps 1 & 2 'Changing the direction of force') and the flexion / extension stops at the rear. Then position the orthotics on the plaster positive and push the locking / unlocking lever to the horizontal position to lock the angle adjustment of the orthotics. Then remove the orthotics from the cast and fit the flexion / extension stops, adjusted to the set angle.

Please note: check that the visible colour on the spring is in accordance with the intended flexion/extension force direction. Place the force output indicator disc with the same colour as the spring visible upwards on the spring cassette, but in such a way that the number 0 is aligned with the proximal force output indicator mark (6) on the brake housing of the system joint. Then insert the washer and housing and secure with the central housing screw. Then replace the spring cassette and release the joint.

Locking the joint

In order to more easily position and remove the orthotics, the joint can be locked at virtually each desired angle. The release and locking mechanism is located to the side of the lower joint half. In order to lock the joint, the catch should be pushed forward/up.

Setting flexion and extension stops

To ensure correct operation, it is **mandatory** to set the flexion and extension stops. Failure to do so may result in excessive strain on the joint. There are two screws at the back of the functional MultiMotion joint to set the range of motion. The standard stop setting in both functional joints is 0° flexion and 95° extension. With the 'Medium' model, the stops can be moved in 15° increments and with the 'Large' model in 12° increments. Use the supplied roundhead screwdriver to loosen the flexion and extension screws, move them as needed and then retighten them.

CAUTION

BEWARE OF RISKS OF INJURY BEFORE UNLOCKING THE JOINT. THE SET SPRING TENSION MUST BE TAKEN INTO ACCOUNT. IF THIS IS HIGH, IT IS ADVISABLE TO HAVE A CARER ASSIST WITH UNLOCKING IN ORDER TO GRADUALLY GUIDE THE ORTHOTICS TOWARDS THE SET FLEXION OR EXTENSION FORCE.

Setting the spring tension

On both sides of the system joint, there is an Allen wrench with which the spring tension can be adjusted.

Please note: There must be no tension on the spring during the fitting phase. It should only be adjusted during rehabilitation, taking into account that the patient may have to get used to wearing the orthotics first. The duration of habituation is determined by the sensibility and is therefore patient-dependent and can vary from hours to days and is determined by the treating specialist/therapist.

Patient documentation form

A patient documentation form is included in the package on which the progression of the active and passive range of motion, spring tension and wearing time can be recorded. Be sure to give this form to the patient/carer or therapist concerned.

Liability

Basko Healthcare as the manufacturer is solely liable if the product is used in accordance with the prescribed conditions and the purpose for which it is intended. Basko Healthcare recommends handling the product in accordance with applicable rules and performing maintenance as described in these user instructions. Non-compliance with the above may cause the MultiMotion joints to not function or function inadequately and may impact the guarantee. Changes, adjustments and modifications to the product not described in these instructions, are the responsibility and liability of the person who performs them.

CE conformity

These products meet the requirements in Directive (EU) 2017/745 for medical devices. These products have been included in Class I based on classification criteria for medical devices as per annex VIII of the Directive. Basko Healthcare therefore wrote up this conformity statement entirely at its own responsibility as per annex XI of the Directive.

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