# Swing Phase Lock 2

Orthotic knee joint

Professional instructions

**Basko**  
Healthcare

[basko.com](http://basko.com)

<b>NL</b>	<b>- Technische gebruiksaanwijzing</b>	<b>3</b>
<b>DE</b>	<b>- Technische Gebrauchsanleitung</b>	<b>15</b>
<b>ENG</b>	<b>- Technical manual</b>	<b>27</b>
<b>FR</b>	<b>- Manuel technique</b>	<b>39</b>
<b>ES</b>	<b>- Instrucciones técnicas</b>	<b>51</b>
<b>IT</b>	<b>- Istruzioni tecniche per l'uso</b>	<b>63</b>
<b>SE</b>	<b>- Tekniska anvisningar</b>	<b>75</b>
<b>DK</b>	<b>- Tekniske instruktioner</b>	<b>87</b>
<b>NO</b>	<b>- Tekniske instruksjoner</b>	<b>99</b>
<b>FI</b>	<b>- Tekniset ohjeet</b>	<b>111</b>

# Nederlands

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# English

User directives	29
Intended use / Indications	29
<b>Indications</b>	29
<b>Contra Indications</b>	29
<b>Functional features / Other terms / Usages</b>	30
<b>Important</b>	30
Function	30
<b>Technical Features SPL 2</b>	30
<b>Technical Features SPC</b>	30
<b>Technical Features Satellite</b>	30
Cleaning / disinfection Instructions	31
Maintenance Instructions	31
Used Materials	31
Constructional guidelines	31
<b>Orthosis construction</b>	31
<b>Measuring</b>	31
<b>Orthosis Construction</b>	32
<b>Joint Alignment</b>	32
<b>Dynamic Initial Fitting</b>	32
<b>Function Check</b>	32
<b>Adjusting locking and unlocking</b>	32
<b>Satellite Control</b>	33
<b>Shortening the Satellite Cable</b>	33
Checking all functions	35
<b>All functions must be checked once again</b>	35
Package contents	35
<b>IQ Swing Phase Lock 2 Set</b>	35
<b>IQ Swing Phase Lock 2 monolateral</b>	35
<b>IQ SPL 2 Joint</b>	35
<b>IQ SPC Joint</b>	35
<b>IQ Satellite</b>	35
Overview Parts	36
<b>Parts SPL 2 Joint</b>	36
<b>Parts SPC Joint</b>	37
<b>Parts Satellite</b>	38

## User directives

**Your doctor is the one whom prescribed the orthosis and stipulates the duration of treatment. The daily use duration prescribed by the doctor must be observed closely. The use of the orthosis is done under the supervision of the doctor. The professional staff, which builds and adapts the articulation system into the orthosis, is responsible for correct functioning of the Joints. In case of skin problems, circulation problems, unusual pain or other complications, which to your opinion are directly related to the product, please contact your doctor or CO at once.**

## Intended use / Indications

Usage of the Swing Phase Lock 2 knee Joint system provides an optimum care for patients with Paresis/ Paralysis (e.g. Polio) and as a consequence contributes to a successful therapeutic treatment. Other usage as prescribed above is not permitted. For a correct application the usage in combination with the satellite in the orthosis is mandatory. The Joint system has been designed only as a non-weight bearing Joint (KAFO/KO) and should not be used in any weight bearing orthoses. Non-observance of these regulations excludes any liability claims. The Joint is available in two different version, for usage as a monolateral Joint or as a bilateral Joint. **When used in a monolateral construction, the maximum patient weight cannot exceed more than 100 kg (220 lbs).**

### Indications

SPL 2 Orthotic Knee Joint bilateral and monolateral:

- |   |   |
|---|---|
| <ul style="list-style-type: none"> <li>• Apoplexy (CVA) (ICD10: I64)</li> <li>• Multiple sclerosis (ICD10: G35)</li> <li>• Myopathy (ICD10: G72.9)</li> </ul> | <ul style="list-style-type: none"> <li>• Paralysis (ICD10: G83.9)</li> <li>• Paresis (peripheral) (ICD10: G83.9)</li> <li>• Poliomyelitis (ICD10: A80)</li> </ul> |
|---|---|

### Contra Indications

SPL 2 Orthotic Knee Joint Bilateral

- Flexion Contracture > 10° (ICD10: M24.59)
- Hip Contracture
- Spasm
- Weight bearing Orthosis

SPL 2 Orthotic Knee Joint Monolateral

- Flexion Contracture > 5°
- Hip Contracture
- Patient weight > 100 kg
- Spasm
- Varus > 10° (ICD10: M21.10/Q74.9)
- Valgus >10° (ICD10: M21.00/Q74.9)
- Weight bearing Orthosis

### Functional features / Other terms / Usages

This articulating system is meant to compensate weakened or Total Failure of the Knee Extension function.

The SPL 2 Joint offers knee stability during stance and offers free movement of the knee during swing. This feature is fully automatically activated with each step. It is prohibited to be used in water or exposed to temperatures exceeding 50°C.

### Important

The patient must be well informed concerning the functions and possibilities of the Knee Joint system.

## Function

### Technical Features SPL 2

The technique of the SPL 2 Joint (Swing Phase Lock) is based on a simple internal pendulum mechanism, which locks and unlocks the knee depending on the angle of the Joint in the sagittal plane. During gait, the device locks at the end of swing phase, just prior to heel strike, and unlocks the knee at heel off in preparation for swing. The Joint cannot be unlocked if it is loaded in flexion. The SPL 2 Joint is always to be mounted on the lateral side and is **suitable for 19 mm Uprights**. The Joint weight is 340 grams.

### Technical Features SPC

De technique of the SPC Joint (Swing Phase Control) usage friction to regulate or influence excessive Knee Flexion. The moment to full flexion is somewhat shortened and slightly variable. The SPC- Joint is always to be mounted on the medial side and is **suitable for 19 mm Uprights**. The Joint weight is 190 grams.


### Technical Features Satellite

The SPL 2 Joint is controlled by a proximal remote push-button switch (Satellite). The Satellite is mounted to the thigh part of the orthosis with the mounting bracket that is supplied in the packaging.

The 3 modes to manually control the SPL 2 Joint are:

1. Automatic locking
2. Permanent unlocking
3. Permanent locking

By deactivating the permanent unlocking mode the Satellite switches into the permanent locking mode and needs to be pushed into the Automatic function.

## Cleaning / disinfection Instructions

The SPL 2 Joint may not be greased, oiled or in any other way lubricated. When cleaned, please use Naphtha (**no acetone, dilutions or turpentine, etc.**) or some similar cleaning fluids or compressed air.

## Maintenance Instructions

Maintenance cycles depend strongly on the patient's activity level. At a moderate activity level, the Joints must be checked / serviced every 6 months. When doing this, attention must be paid to wear and tear, Joint play and damages. If the Joint is no longer running freely, a thorough service is recommended.

If the Joint has been, due to circumstances, extremely heavily loaded in a flexed position e.g. with stumbling, tripping or walking (falling) down the staircase, the Joint has to be inspected by the Orthotist. In case the patient constantly walks with a slightly flexed knee or regularly usage the safety stops, adjustment of the extension position towards flexion is absolutely necessary.

## Used Materials

SPL 2 Joint: Housing: Stainless Steel; Function Unit: brass; Cover plates: ABS plastic

SPC Joint: Housing: Stainless Steel; Cover plates and Bumper Pin: ABS plastic

Satellite: Housing: ABS, POM plastic; Cable: Teflon

## Constructional guidelines

### Orthosis construction

#### Measuring

Casting must be done with a fully, in neutral, extended leg. Stretching however should be possible without too much effort to allow full extension (needed to unlock the Joint) once the Orthosis is worn. When the Alignment Unit is used (**Art. nr. AK-100/AK-101**) to construct the Orthosis, make sure the square with the T-bar is placed horizontal in the Frontal Plane to keep the axis straight. Alignment should be based on the plumb line to create a neutral baseline for the Function Unit.

## Orthosis Construction

Basically any type of construction is suitable as long as sufficient **torsion rigidity** and a **possibility of changing the flexion position** during the initial fitting is taken into account.

When the Joint is used in a Carbon fiber reinforced Construction, the SPL 2 should be disassembled by removing the Satellite – disconnect the cable within the Function Unit (**H**) and loosen the IQ set screws (**D**), Loosen the adjustment screw and gently pull the cable from the housing. Now the Function Unit (**H**) can be pushed inside and removed.

The SPC flexion Control Unit has to be removed (**E**) and sealed with Stick Wax (**Art. nr. 990035**) to prevent improper resin flow. In this case the Upright spacer can be used (**Art. nr. 029860**) and aligned, allowing the Joints to be mounted in a later stage.

## Joint Alignment


The Alignment Unit has to be parallel to the plumb line. The SPC Joint should not be aligned in full extension. If the Joint is aligned without using the Unit the Cover Plates (**B**) may be removed allowing the Joint axis to be used for parallel alignment.

## Dynamic Initial Fitting

### Function Check

Prior to the initial fitting the proper setting of the Joint need to be checked

- The slotted arm of the Function Unit (**H**) must be set in neutral 0 position (**see picture**) when the Satellite is in the 1st mode. Adjust with the Proximal adjustment screw (**F**) on the SPL 2 is necessary.
- Can the SPL 2 Joint lock after swing with the Satellite in the 1st mode?

### Adjusting locking and unlocking

Timing settings to lock or unlock can be read or adjusted by moving the pendulum weight of the Function Unit.

- Moving the pendulum ANTERIOR (forward), DECREASES the tendency to lock, because its weight is shifted forward.
- Moving the pendulum POSTERIOR (backward), INCREASES the tendency to lock because its weight is shifted backward.

## Satellite Control

The patient needs to be instructed on functionality and how to operate the satellite.

Fig. 1

### Mode 1 – Automatic Locking and Unlocking

The slider on the satellite is positioned in the middle. This is the normal position by which the SPL 2 joint automatically locks and unlocks during walking. In this mode, the joint can also be unlocked once, by slightly pushing the switch upwards. The automatic mode is activated again when the joint is fully extended.

Fig. 2

### Mode 2 – Permanently Unlocked

The slider on the satellite is slid upward from the center. In this mode the SPL 2 joint is permanently unlocked. **NOTE:** When deactivating the permanent unlocking (by pushing the central button), the satellite automatically shifts into the permanent locking mode and thus needs to be slid up into the automatic locking/unlocking mode.

Fig. 3

### Mode 3 – permanently Locked

For this purpose the central button on the satellite is to be pressed for the slider to slid fully down. This is a safety mode in which the SPL 2 joint is constantly locked (e.g. at home, making many turns, etc.).

## Shortening the Satellite Cable

First determine the position where the satellite will be mounted on the thigh cuff of the orthosis in order to establish the correct cable length. Bear in mind that the cable should have a slight S-shape.

1. Position the joint in full extension and set the satellite in Mode 3 by pressing the button on the satellite (see Fig. 3). Make sure the joint is indeed locked.
2. Use the supplied Allen screwdriver to loosen the two screws on the medial (stainless steel) Cover Plate and remove it (see Fig. 4).
3. After removal of the cover plate check if the upper side of the Lever Arm (holding the satellite cable) is parallel to the lower of the three markers on the Function Unit (see Fig. 5).
4. Using a Phillips screwdriver, loosen the screw of the Service Hatch at the back of the satellite and remove the cover and screw.
5. Use the Allen screwdriver to loosen the cone point safety screw two to three turns and slide the Cable End Block of the inner cable (see Fig. 6).
6. Twist the Teflon outer cable from the Swivel Connector at the bottom of the satellite. Place the Swivel in one hand between your thumb and forefinger and the Teflon outer cable in the other hand between the thumb and forefinger and unscrew (counterclockwise) the outer cable from the Swivel Connector (see Fig. 7). Subtract the inner cable outwards through the Swivel.
7. Loosen the Cone Point Safety Screw on the side of the joint, two turns will do (see Fig. 8). Do not remove the screw from the joint.
8. Turn the Satellite Set Screw (which connects the Teflon cable to the joint) completely in the joint and then two complete turns up again.
9. **NOTE:** carve the Teflon outer cable with a sharp knife to the desired length and remove the excess outer cable (the inner steel cable will be shortened in step 14).


Fig. 4


Fig. 5


Fig. 6

10. Now slide the inner cable through the bottom of the Swivel Connector into the Satellite Housing and screw the Teflon outer cable onto the Swivel. Place the Swivel in one hand between your thumb and forefinger and the Teflon cable in the other hand between the thumb and forefinger. Press and turn (clockwise) the two parts firmly together (the Swivel is self-tapping).
11. Slide the Cable End Block back over the steel inner cable and position the End Block back into the corresponding slot in the satellite (see Fig. 9).
12. Recheck the correct position (Mode 3) of the Lever Arm in the Function Unit of the joint (see Fig. 5).
13. Tighten the Cone Point Safety screw of the Cable End Block firmly on the inner cable.
14. Use a sharp cutter to cut off the excess inner cable at  $\pm 5$  mm above the Cable End Block.
15. Install the Service Hatch back on the back of the satellite and reattach the screw.
16. Slide the satellite switch into the Automatic Position (Mode 1) (see Fig. 1) and check within the Function Unit of the joint if the top of the Lever Arm has moved and now is parallel with the center marker (see Fig. 10).
17. Use the Satellite Set Screw (which connects the Teflon cable to the joint) to fine tune and adjust the position described in step 16.
18. Tighten the Cone Point Safety screw back in the joint (see Fig. 8).
19. Mount the medial Cover Plate back onto the joint and tighten the screws.
20. Loosen both screws on the front side of the Satellite Housing and remove the Auxiliary Bracket.
21. Determine the exact position where the satellite is to be mounted on the thigh cuff and contour the bracket if necessary. Mount the Auxiliary Bracket on the brace and reattach the satellite using the screws at the front again.
22. Use the provided Cable Clamps to attach the satellite cable in order to ensure an optimal, flat resting the cable run against the thigh cuff of the orthosis.


Fig. 7


Fig. 8


Fig. 9


Fig. 10

Once all the adjustments are done, the opening on the Function Unit must be re-sealed using the IQ sticker. This reduces the likelihood of dirt and clothing fibers getting in.

## Checking all functions

### All functions must be checked once again

- Do all three modes on the satellite switch reliably?
- Can a reliable locking be achieved with free swing?
- Is the locking sufficient (is the cable to the satellite on the right length, is the inner cable in a neutral position, in the 0 position)?
- Does the Joint unlock reliably after the back swing?
- Are the stops all again in position and does the Joint run without play and smoothly?

## Package contents

### **IQ Swing Phase Lock 2 Set**

(Art.nr.: IQ200R or L)

- 1 SPL 2 Joint
- 1 SPC Joint
- 1 Satellite and cable
- 1 Inner hexagonal wrench
- 2 Stickers
- 2 Cable Clamps
- 2 Nyloplex stud with header
- 1 Manual

### **IQ SPL 2 Joint**

(Art.nr.: IQ210R or L)

- 1 SPL 2 Joint
- 1 Inner hexagonal wrench
- 2 Stickers
- 2 Cable Clamps
- 2 Nyloplex stud with header
- 1 Manual

### **IQ Satellite**

(Art.nr.: IQ141)

- 1 Satellite and Cable
- 1 Manual

### **IQ Swing Phase Lock 2 monolateral**

(Art.nr.: IQ250R or L)

- 1 SPL 2 Joint
- 1 Satellite and cable
- 1 Inner hexagonal wrench
- 2 Stickers
- 2 Cable Clamps
- 2 Nyloplex stud with header
- 1 Manual

### **IQ SPC Joint**

(Art.nr.: IQ120R or L)

- 1 SPC Joint
- 1 Inner hexagonal wrench
- 1 Manual

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