

# Carbon Footplate Rigid

Hiilikuitupohjallisen tarkoituksena on suojata metatarsophalangeaalinivelet MTP-nivelet ja jalkaterän etuosaa rajoittamalla liikelaajuutta (ROM). Hiilikuitupohjallista voidaan käyttää myös immobilisaatiohoidon jälkeen POST-OP-kengässä tai diabetesjalkineessa liikkumisen asteittaiseksi lisäämiseksi ja tasapainon parantamiseksi.

Hiilikuitupohjallinen auttaa jakamaan jalkaterän etuosaan kohdistuvaa painetta ja helpottaa ylipronaatiota. Hiilikuitupohjallinen on peitettävä kengän omalla tai yksilöllisesti valmistetulla pohjallisella.



## Käyttöalueet

Hallux rigidus, hallux limitus, rasitusmurtumat, etujalkaterän vammat, varpaiden nyrjähdykset ja muut etujalkaterän/keskijalkaterän vammat.

## Kontraindikaatiot

Avoimet haavaumat jalassa, vakavat murtumat ja yli 135 kg:n paino.

Tuotenro.	Malli	Paksuus	Koko	Vas/Oik	Pituus	Leveys
280361011	Carbon Footplate Rigid	2,0 mm	Small	Vasen	235 mm	77 mm
280361012	Carbon Footplate Rigid	2,0 mm	Medium	Vasen	255 mm	80 mm
280361013	Carbon Footplate Rigid	2,0 mm	Large	Vasen	275 mm	90 mm
280361014	Carbon Footplate Rigid	2,0 mm	X-Large	Vasen	305 mm	99 mm
280362011	Carbon Footplate Rigid	2,0 mm	Small	Oikea	235 mm	77 mm
280362012	Carbon Footplate Rigid	2,0 mm	Medium	Oikea	255 mm	80 mm
280362013	Carbon Footplate Rigid	2,0 mm	Large	Oikea	275 mm	90 mm
280362014	Carbon Footplate Rigid	2,0 mm	X-Large	Oikea	305 mm	99 mm